


COMPACT Automatic Soft Tissue Biopsy Device

One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer.

- Programmable: releases the stylet allowing time for tissue to settle in notch (ideal for denser tissue or when there is a desire to verify notch position) and activates the cutting cannula
- Thin-wall cannula for larger core
- Colour coded for quick gauge size identification
- Depth markings etched on cannula for precise positioning
- The Co-Axial kit permits multiple biopsies (US and CT guided)
- Loaded status indicator
- Plastic safety protection to avoid accidental firing



Co-Axial kit

COMPACT

PRODUCT CODE:	SIZE	CO-AXIAL KIT
● VCT1210	12G x 10cm	VCTC 1210
● VCT1410	14G x 10cm	VCTC 1410
● VCT1411	14G x 11cm	VCTC 1411
● VCT1413	14G x 13cm	VCTC 1413
● VCT1415	14G x 15cm	VCTC 1415
○ VCT1610	16G x 10cm	VCTC 1610
○ VCT1615	16G x 15cm	VCTC 1615
○ VCT1620	16G x 20cm	VCTC 1620
● VCT1810	18G x 10cm	VCTC 1810
● VCT1815	18G x 15cm	VCTC 1815
● VCT1820	18G x 20cm	VCTC 1820
● VCT1825	18G x 25cm	VCTC 1825
● VCT2010	20G x 10cm	VCTC 2010
● VCT2015	20G x 15cm	VCTC 2015

(Box of x10 pieces)



One handed design with loaded status indicator



Colour code

Needles with different diameters and lengths may be provided on request.