

COREX Automatic Biopsy Device

- One-handed design simplifies loading and operation whilst allowing for convenient manipulation of the ultrasound transducer
- Equipped with a visual and mechanical safety-system
- Sample notch of 19 mm
- Thin-wall cannula for larger core
- Two triggers permit ergonomic handling for right and left-handed application, regardless of the position of the extraction location
- Colour coded for quick gauge size identification
- Depth markings etched on cannula for precise positioning
- The Coaxial kit permits multiple biopsies (US and CT guided)



COREX

PRODUCT CODE:	SIZE	CO-AXIAL KIT
● VBS 1410	14G x 10cm	VBSC 1410
● VBS 1415	14G x 15cm	VBSC 1415
○ VBS 1610	16G x 10cm	VBSC 1610
○ VBS 1615	16G x 15cm	VBSC 1615
○ VBS 1620	16G x 20cm	VBSC 1620
○ VBS 1625	16G x 25cm	VBSC 1625
● VBS 1810	18G x 10cm	VBSC 1810
● VBS 1815	18G x 15cm	VBSC 1815
● VBS 1820	18G x 20cm	VBSC 1820
● VBS 1825	18G x 25cm	VBSC 1825

(Box of x10 pieces)

Needles with different diameters and lengths may be provided on request.

Mechanical safety-system



Longer handle for ergonomic improvement

