

## INTRO/INTRO PLUS Standard C-Axial Introducer Needle used for Single Access in multiple Biopsies

- Echogenic treatment for ultrasound guided procedures
- Centimetre marking and optional sliding stopper
- Sterilised by ETO, shelf life 5 years
- INTRO (code VVT) with blunt cannula and trocar tip
- INTRO PLUS (code WTC) with sharp cutting cannula and trocar tip



### INTRO/INTRO PLUS

| PRODUCT CODES:                                   |           | SIZE       |
|--|-----------|------------|
| <span style="color: green;">●</span> VVT1105     | VVTC 1105 | 11G x 5cm  |
| <span style="color: lightblue;">●</span> VVT1305 | VVTC 1305 | 13G x 5cm  |
| <span style="color: lightblue;">●</span> VVT1310 | VVTC 1310 | 13G x 10cm |
| <span style="color: lightblue;">●</span> VVT1315 | VVTC 1315 | 13G x 15cm |
| <span style="color: blue;">●</span> VVT1505      | VVTC 1505 | 15G x 5cm  |
| <span style="color: blue;">●</span> VVT1510      | VVTC 1510 | 15G x 10cm |
| <span style="color: blue;">●</span> VVT1515      | VVTC 1515 | 15G x 10cm |
| <span style="color: red;">●</span> VVT1705       | VVTC 1705 | 17G x 5cm  |
| <span style="color: red;">●</span> VVT1710       | VVTC 1710 | 17G x 10cm |
| <span style="color: red;">●</span> VVT1715       | VVTC 1715 | 17G x 15cm |

(Box of x20 pieces)

Depth marker with code colour



Needles with different diameters and lengths may be provided on request.