

RELOCK PREMIUM SOFT Preoperative Repositionable Breast Localisation Device

One-handed deployment for convenient manipulation of the ultrasound transducer.

- Preoperative marking of local lesions
- Retractable hook with unique flexibility of its polyamide guiding thread
- Indicator on wing shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimetre depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- 9mm loop
- Single packaged stiffener cannula used in the operating room




Full Exposure



Half Exposure



Stiffener Cannula



RELOCK PREMIUM SOFT

PRODUCT CODE:	SIZE
● VVPRS1905	19G x 5cm
● VVPRS1907	19G x 7cm
● VVPRS1910	19G x 10cm
● VVPRS2005	20G x 5cm
● VVPRS2007	20G x 7cm
● VVPRS2010	20G x 10cm

(Box of x10 pieces)



Mammography

Needles with different diameters and lengths may be provided on request.