

SURGIMARC SLIM One-Handed Preoperative Breast Localisation Device


- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Puncture cannula with ultrasonic tip with guidance wire in place
- The Z-shaped hook provides safer delivery in critical areas (available in 18G version)
- The X-shaped hook provides excellent fixation and is intended to limit migration (available in 20 and 21G version)



Single Hook



Z-Shaped Hook



X-Shaped Hook



One-Handed Design

SURGIMARC SLIM

PRODUCT CODES:			SIZE:
SINGLE HOOK	X-SHAPED DOUBLE ANCHORING HOOK	Z-SHAPED HOOK	
● VVNS1807		VVNSX1807	18G x 7cm
● VVNS1810		VVNSX1810	18G x 10cm
● VVNS2007	VVNSX2007		20G x 7cm
● VVNS2010	VVNSX2010		20G x 10cm
● VVNS2015	VVNSX2015		20G x 15cm
● VVNS2107	VVNSX2107		21G x 07cm
● VVNS2110	VVNSX2110		21G x 10cm

(Box of x20 pieces)

Needles with different diameters and lengths may be provided on request.