

SURGIMARC Preoperative Breast Localisation Device


- Preoperative marking of local lesions
- Puncture cannula with ultrasonic tip with guidance wire in place
- Transparent cannula connector
- The Z-shaped hook provides safer delivery in critical areas
- The X-shaped hook provides excellent fixation and is intended to limit migration



Single Hook



Z-Shaped Hook



X-Shaped Hook



SURGIMARC

PRODUCT CODES:			SIZE:
● VVN 1807	VVNX1807	VVNZ1807	18G x 7cm
● VVN 1810	VVNX1810	VVNZ1810	18G x 10cm
● VVN 1815	VVNX1815	VVNZ1815	18G x 15cm
● VVN 2005	VVNX2005	VVNZ2005	20G x 5cm
● VVN 2007	VVNX2007	VVNZ2007	20G x 7cm
● VVN 2009	VVNX2009	VVNZ2009	20G x 9cm
● VVN 2010	VVNX2010	VVNZ2010	20G x 10cm
● VVN 2015	VVNX2015	VVNZ2015	20G x 15cm
● VVN 2107	VVNX2107	VVNZ2107	21G x 07cm
● VVN 2110	VVNX2110	VVNZ2110	21G x 10cm

(Box of x20 pieces)

Needles with different diameters and lengths may be provided on request.