

CYTOCAN Chiba Needle for Aspiration Biopsy

- Puncture cannula developed especially for ultrasonic diagnostics
- Echogenic needle tip assists precise needle tip positioning
- Safe for patient and user, thanks to the adjustable sliding depth stopper
- Graduated puncture cannula to assist in precise depth placement of the needle
- Transparent luer cone for visual inspection of aspirate
- Colour coded hub for easy identification of needle gauge size



CYTOCAN

PRODUCT CODES:	SIZE
● VVB 1810	18G x 10cm
● VVB 1815	18G x 15cm
● VVB 1820	18G x 20cm
● VVB 2010	20G x 10cm
● VVB 2015	20G x 15cm
● VVB 2020	20G x 20cm
● VVB 2105	21G x 5cm
● VVB 2110	21G x 10cm
● VVB 2115	21G x 15cm
● VVB 2120	21G x 20cm
● VVB 2210	22G x 10cm
● VVB 2215	22G x 15cm
● VVB 2220	22G x 20cm
● VVB 2310	23G x 10cm
● VVB 2315	23G x 15cm
● VVB 2505	25G x 5cm

(Box of x20 pieces)

Needles with different diameters and lengths may be provided on request.

Reusable syringe gun (code: SG01)

