

## MYELO Bone and Bone Marrow Biopsy Device

- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely



### MYELO

#### VVIE<sub>xyy/xx(x)</sub> for bone marrow harvesting needle

PRODUCT CODES:		SIZE:
VVI11410/50	VVIE1410/50	14G x 30cm
VVI11430/70	VVIE1430/70	14G x 50cm
VVI11510/50	VVIE1510/50	15G x 30cm
VVI11530/70	VVIE1530/70	15G x 50cm
VVI11610/50	VVIE1610/50	16G x 30cm
VVI11630/70	VVIE1630/70	16G x 50cm
VVI11810/50	VVIE1810/50	18G x 30cm

(Box of x10 pieces)

Needles with different diameters and lengths may be provided on request.

Side distal holes VVIEI<sub>xyy/xx(x)</sub>

