


## ERGO Bone and Bone Marrow Biopsy Device

- Ergonomic T-handle design to minimize hand and arm fatigue
- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely



### MYELO


#### VVILExxy/xx(x) for bone marrow harvesting needle

PRODUCT CODES:		SIZE:
VVIL1410/50	VVILE1410/50	14G x 30cm
VVIL1430/70	VVILE1430/70	14G x 50cm
VVIL1510/50	VVILE1510/50	15G x 30cm
VVIL1530/70	VVILE1530/70	15G x 50cm
VVI1610/50	VVILE1610/50	16G x 30cm
VVI1810/50	VVILE1810/50	18G x 30cm

(Box of x10 pieces)

Needles with different diameters and lengths may be provided on request.

Side distal holes VVILExxy/xx(x)



Colour code

