

MYELO Bone and Bone Marrow Biopsy Device

- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely



MYELO

VVIExxyy/xx(x) for bone marrow harvesting needle		
PRODUCT CODES:		SIZE:
VVI1410/50	VVIE1410/50	14G x 30cm
VVI1430/70	VVIE1430/70	14G x 50cm
VVI1510/50	VVIE1510/50	15G x 30cm
VVI1530/70	VVIE1530/70	15G x 50cm
VVI1610/50	VVIE1610/50	16G x 30cm
VVI1630/70	VVIE1630/70	16G x 50cm
VVI1810/50	VVIE1810/50	18G x 30cm

(Box of x10 pieces)

Needles with different diameters and lengths may be provided on request.

Side distal holes VVIExxyy/xx(x)

