

B CUT Bone and Bone Marrow Biopsy Device

- Designed for quick and easy cortical penetration
- Ergonomic T-handle for increased control
- Stylet with pyramidal tip
- Cannula with fish mouth tip or with crown tip
- Centimetre depth markings assist in needle placement
- Improved patient comfort
- Extraction cannula with centimetre markings and blunt tip



B-CUT

VVJxxy: Fish mouth tip collection cannula				
VVJxxy10: Side distal hole for bone marrow aspiration				
VVJxxy20: Crown tip cannula				
VVJExxy: Side distal holes for efficient aspiration during harvesting procedures				
PRODUCT CODES:				SIZE:
● VVJ0810	VVJ081010	VVJ081020		8G x 10cm
● VVJ0815	VVJ081510	VVJ081520		8G x 15cm
● VVJ0910	VVJ091010	VVJ091020		9G x 10cm
● VVJ0915	VVJ091510	VVJ091520		9G x 15cm
○ VVJ1110	VVJ111010	VVJ111020	VVJE1110	11G x 10cm
○ VVJ1115	VVJ111510	VVJ111520	VVJE1115	11G x 15cm
● VVJ1307	VVJ130710	VVJ130720		13G x 7cm
● VVJ1310	VVJ131010	VVJ131020	VVJE1310	13G x 10cm

(Box of x10 pieces)

Side distal holes (VVJExxy)



Fish mouth tip

