


HYSTO SYSTEM II Disposable Automatic System

With two programmable firing modes, these lightweight automated systems give you the precise control and quality sampling capability you need when working with calcified or fibrous lesions. Available with or without coaxial introducer needle. It can be easily operated with one hand, making it ideal for procedures requiring ultrasound guidance

- Large 20 mm sample notch
 - Echogenic tip for accurate placement under ultrasound guidance
 - Two fire trigger and one end operation
 - Can be used with echogenic Coaxial introducer needles
 - One-handed operation for convenient manipulation of guidance device (ultrasound transducer)
- Manipulation of guidance device (ultrasound transducer)



HYSTO SYSTEM II Coaxial Biopsy System

Allows multiple sampling from one percutaneous puncture with minimal damage to surrounding tissue.

- Trocar point tip maintains directional stability
- Blunt tip cannula helps prevent inadvertent distal drift and may reduce risk of tissue trauma
- Echogenic markings facilitate precise ultrasound positioning
- Locking skin stop helps prevent forward migration
- Matched in length to cutting cannula



HYSTO SYSTEM II Disposable Automatic Biopsy Needle

PRODUCT CODE:	DESCRIPTION
HSB 1410	14G x 10cm
HSB 1416	14G x 16cm
HSB 1420	14G x 20cm
HSB 1610	16G x 10cm
HSB 1616	16G x 16cm
HSB 1620	16G x 20cm
HSB 1810	18G x 10cm
HSB 1816	18G x 16cm
HSB 1820	18G x 20cm
HSB 1825	18G x 25cm
HSB 2010	20G x 10cm
HSB 2016	20G x 16cm
HSB 2020	20G x 20cm

(Box of x10 pieces)

Coaxial Biopsy Needle for Hysto System II

PRODUCT CODE:	DESCRIPTION
CPX 1410	13G x 5cm Coaxial to HSB 1410
CPX 1416	13G x 11cm Coaxial to HSB 1416
CPX 1420	13G x 15cm Coaxial to HSB 1420
CPX 1610	15G x 5cm Coaxial to HSB 1610
CPX 1616	15G x 11cm Coaxial to HSB 1616
CPX 1620	15G x 15cm Coaxial to HSB 1620
CPX 1810	17G x 5cm Coaxial to HSB 1810
CPX 1816	17G x 11cm Coaxial to HSB 1816
CPX 1820	17G x 15cm Coaxial to HSB 1820
CPX 2010	19G x 5cm Coaxial to HSB 2010
CPX 2016	19G x 11cm Coaxial to HSB 2016
CPX 2020	19G x 15cm Coaxial to HSB 2020

(Box of x10 pieces)