

Fine Needle Aspiration (FNA) Biopsy Needle

Fine needle aspiration for cytological and histological evaluation.

- Sharpened cannula for a safe and painless insertion and penetration
- Available with echogenic tip, ensuring accurate ultrasound placement
- Sliding stop and centimetre division promoting accurate depth placement
- Centimetre marking and color coded for identification



Chiba Biopsy Needle

PRODUCT CODE:	DESCRIPTION
CNA 1809	18G x 9cm
CNA 1815	18G x 15cm
CNA 1820	19G x 20cm
CNA 2009	20G x 9cm
CNA 2015	20G x 15cm
CNA 2020	20G x 20cm
CNA 2209	22G x 9cm
CNA 2215	22G x 15cm
CNA 2220	22G x 20cm
CNA 2309	23G x 9cm
CNA 2509	25G x 9cm
CNA 2009E	20G x 9cm
CNA 2209E	22G x 9cm
CNA 2020E	20G x 20cm
CNA 2220E	22G x 20cm

(Box of x20 pieces)

E Codes are without cm markings or sliding stopper



Franseen Biopsy Needle

PRODUCT CODE:	DESCRIPTION
FRA 1809	18G x 9cm
FRA 1815	18G x 15cm
FRA 1820	19G x 20cm
FRA 2009	20G x 9cm
FRA 2015	20G x 15cm
FRA 2020	20G x 20cm
FRA 2209	22G x 9cm
FRA 2215	22G x 15cm
FRA 2220	22G x 20cm

(Box of x10 pieces)



Westcott Biopsy Needle

PRODUCT CODE:	DESCRIPTION
WES 1809	18G x 9cm
WES 1815	18G x 15cm
WES 1820	19G x 20cm
WES 2009	20G x 9cm
WES 2015	20G x 15cm
WES 2020	20G x 20cm
WES 2209	22G x 9cm
WES 2215	22G x 15cm
WES 2220	22G x 20cm

(Box of x10 pieces)



Spinal Needle

PRODUCT CODE:	DESCRIPTION
SNA 1809	18G x 9cm
SNA 1909	19G x 9cm
SNA 2009	20G x 9cm
SNA 2109	21G x 9cm

(Box of x10 pieces)