

RETRIEVE Jamshidi Bone Marrow Biopsy System


RETRIEVE Bone Marrow Biopsy System is designed to provide quality biopsies following the classical Jamshidi procedure. It is also available in the explant/transplant version.

- Cannula with 5 cutting edges to facilitate penetration
- Stylet with trocar point for smoother entry into the marrow
- Luer lock hub
- Provided with a special cap for a safe removal of the sample
- Available with Explant Version: 5 lateral holes on the external cannula and additional stylet (Code RDSE)

RETRIEVE Bone Marrow Biopsy

PRODUCT CODE:	DESCRIPTION
RDS 0810	8G x 10cm
RDS 0815	8G x 15cm
RDS 0910	9G x 10cm
RDS 0915	9G x 15cm
RDS 1110	11G x 10cm
RDS 1115	11G x 15cm
RDS 1307	13G x 7cm
RDS 1310	13G x 10cm
RDSE 1110	11G x 10cm, Explant/transplant version
RDSE 1310	13G x 10cm, Explant/transplant version

(Box of x5 pieces)



Code RDS

Code RDSE Tip with 5 side holes.