

## ISTAN Bone Marrow Aspiration System

ISTAN Bone Marrow Aspiration System to be used either in the sternum or iliac crest. It is also available in the explant/transplant version.

- Triple sharpened tip cannula to allow a non traumatic, fast and easy penetration into the bone marrow cavity
- An adjustable, removable needle stop safely controls depth of penetration
- Luer lock hub
- Available with Explant/Transplant Version: 2 lateral holes on the external cannula (Code ISTE)

### ISTAN Bone Marrow Aspiration

PRODUCT CODE:	DESCRIPTION
IST1501	15G x 10–44mm
IST1503	15G x 25–60mm
IST1601	16G x 10–44mm
IST1603	16G x 25–60mm
IST1801	18G x 10–44mm
IST1803	18G x 25–60mm
ISTE1404	14G x 37–71mm, Explant/transplant version

(Box of x10 pieces)



Code IST



Code ISTE with  
2 side holes