

## BMH® Bone Marrow Aspiration Needle

BMH® Bone marrow Aspiration system to be used either in the sternum or iliac crest. Ergonomic handle allow a better grip during procedure.

- Triple sharpened tip cannula
- Luer lock hub
- Ergonomic handle for better grip during procedure
- Available with Explant/Transplant Version: 2 lateral holes on the external cannula (Code BMHE)

### BMH Bone Marrow Aspiration

PRODUCT CODE:	DESCRIPTION
BMH1110	11G x 10cm
BMH1115	11G x 15cm
BMH1310	13G x 10cmm
BMH1507	15G x 7cm
BMH1510	15G x 10cm
BMH1512	15G x 12cm
BMHE1507	15G x 7cm, Explant/transplant version

(Box of x10 pieces)



Code BMH



Code BMHE  
with 2 side  
holes