

## ISG Hysterosalpingography Catheter

ISG facilitates injection of contrast media into the uterine cavity to evaluate patency of the fallopian tubes. The catheter design allows good physician control while maintaining sufficient flexibility for comfort.

- The rounded catheter tip and large side holes assures safe placement and superior flow
- Choice of two sizes 5Fr and 7Fr
- Balloon latex free
- The catheter is supplied with introducer, stopcock and syringe



### ISG Hysterosalpingography Catheter

PRODUCT CODE:	DESCRIPTION
76005	5Fr x 30cm
76007	7Fr x 30cm

(Box of x5 pieces)

