

HERNIATOME


Herniectomy with Herniatome, performed with CT and Fluoroscopy or fluoroscopy alone is a minimally invasive procedure on the spine, the aim of which is identical to that of surgery to extract the hernia or a portion of it to reduce the pressure on the nerve root.

- Evaluation in real time of the decompression
- Qualitative and quantitative evaluation
- Cannula with lateral window, that allow a better decompression and extraction of the hernia by covering a wider area with minimal movement
- Available with straight and curved Tip
- Available with 20ga introducer needle that ensures the safe and correct positioning of the Herniatome

HERNIATOME - Curved Tip (Direct Access Approach)

PRODUCT CODE:	DESCRIPTION
5091709	17G x 9cm, Cervical Type with 20ga - 15cm Introducer
5091715	17G x 15cm, Lumber Type with 20ga - 30cm Introducer
5091809	18G x 9cm, Cervical Type with 20ga - 15cm Introducer

(Box of x1 piece)



HERNIATOME - Straight Tip (Direct Access Approach)

PRODUCT CODE:	DESCRIPTION
5091749	17G x 9cm, Cervical Type
5091745	17G x 15cm, Lumber Type
5091849	18G x 9cm, Cervical Type

(Box of x1 piece)

