

P-DRAIN

Pigtail Universal Drainage Catheter.

- Radiopaque polyurethane catheter
- Stiffened cannula with trocar stylet
- Available with stringlock and hydrophilic coating
- Large oval holes in pigtail for maximum drainage
- Position marking reference catheter migration post-placement
- Hydrophilic coating to reduce friction during insertion



'P-Drain Hydrophilic Universal Drainage Catheter with STRINGLOCK'



P-DRAIN Universal Drainage Catheter

PRODUCT CODE:	DESCRIPTION
14/0630	6Fr x 30cm
14/0730	7Fr x 30cm
14/0830	8Fr x 30cm
14/0930	9Fr x 30cm
14/1030	10Fr x 30cm
14/1230	12Fr x 30cm

(Box of x5 pieces)

P-DRAIN Universal Drainage Catheter SINGLE REPLACEMENT Catheter

PRODUCT CODE:	DESCRIPTION
15/0630	6Fr x 30cm
15/0730	7Fr x 30cm
15/0830	8Fr x 30cm
15/0930	9Fr x 30cm
15/1030	10Fr x 30cm
15/1230	12Fr x 30cm

(Box of x5 pieces)

P-DRAIN Hydrophilic Universal Drainage Catheter with STRINGLOCK

PRODUCT CODE:	DESCRIPTION
34/0630	6Fr x 30cm
34/0730	7Fr x 30cm
34/0830	8Fr x 30cm
34/0930	9Fr x 30cm
34/1030	10Fr x 30cm
34/1230	12Fr x 30cm

(Box of x5 pieces)

P-Drain Hydrophilic Universal Drainage Catheter SINGLE REPLACEMENT Catheter

PRODUCT CODE:	DESCRIPTION
25/0630	6Fr x 30cm
25/0730	7Fr x 30cm
25/0830	8Fr x 30cm
25/0930	9Fr x 30cm
25/1030	10Fr x 30cm
25/1230	12Fr x 30cm

(Box of x5 pieces)

P-Drain Hydrophilic Universal Drainage Catheter

PRODUCT CODE:	DESCRIPTION
24/0630	6Fr x 30cm
24/0730	7Fr x 30cm
24/0830	8Fr x 30cm
24/0930	9Fr x 30cm
24/1030	10Fr x 30cm
24/1230	12Fr x 30cm

(Box of x5 pieces)