

B LOCK Bone and Bone Marrow Biopsy Device

- Designed for quick and easy cortical penetration
- Ergonomic T-handle for increased control
- Stylet with pyramidal tip
- Cannula with fish mouth tip or with crown tip
- Centimetre depth markings assist in needle placement
- Spoon collection cannula or forceps collection cannula
- Improved patient comfort



B-LOCK

VVJTxyy: Fish mouth tip and Spoon collection cannula		
VVJTxyy10: Crown tip and Forceps collection cannula		
PRODUCT CODE:		SIZE:
● VVJT0810	VVJT081010	8G x 10cm
● VVJT0815	VVJT081510	8G x 15cm
● VVJT0910	VVJT091010	9G x 10cm
● VVJT0915	VVJT091510	9G x 15cm
○ VVJT1110	VVJT111010	11G x 10cm
○ VVJT1115	VVJT111510	11G x 15cm
● VVJT1307	VVJT130710	13G x 7cm
● VVJT1310	VVJT131010	13G x 10cm

(Box of x10 pieces)

Needles with different diameters and lengths may be provided on request.

Spoon collection cannula



Forceps collection cannula

