

MYELO SLIM Bone and Bone Marrow Biopsy Device

- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely



MYELO SLIM

VVIE1xxyy/xx(x) for bone marrow harvesting needle		
PRODUCT CODES:		SIZE:
VVI11410/50	VVIE11410/50	14G x 30cm
VVI11430/70	VVIE11430/70	14G x 50cm
VVI11510/50	VVIE11510/50	15G x 30cm
VVI11530/70	VVIE11530/70	15G x 50cm
VVI11610/50	VVIE11610/50	16G x 30cm
VVI11630/70	VVIE11630/70	16G x 50cm
VVI11810/50	VVIE11810/50	18G x 30cm

(Box of x10 pieces)

Needles with different diameters and lengths may be provided on request.

Side distal holes VVIE1xxyy/xx(x)

