

# COVER-ENT®

## STERILE SINGLE-USE SHEATHS FOR FLEXIBLE NASOPHARYNGOSCOPES



Made  
in Italy



Quality Control  
on 100% of  
production



+ Maintains the quality of  
vision

Sterile, practical, and  
ready to use

Greater safety and  
productivity

Effective barrier against  
microorganisms

Lower risks of  
infection and  
adverse events

Preserving the efficiency  
of nasopharyngoscopes

Reducing maintenance  
costs and times

# COVER

CE 0477

Regulation (EU) 2017/745

Cover S.r.l, founded in 2005, was born from the experience of entrepreneurs working in specific sectors since the 80's with the mission to develop and produce high quality medical devices. The goal has been achieved thanks to our long experience and following our clients' suggestions.

# COVER-ENT®

## STERILE SINGLE-USE SHEATHS FOR FLEXIBLE NASOPHARYNGOSCOPES



### PACKAGING:

20 pcs per package,  
5 packs per box

WE ARE ABLE TO PRODUCE SHEATHS FOR ANY  
TYPE OF NASOPHARYNGOSCOPE NOT IN THE  
BROCHURE



Certified Quality Management System  
UNI CEI EN ISO  
13485:2021



Distributed by:

Manufactured by:

# COVER

A YUSHU Healthcare Company

YUSHU  
HEALTHCARE

ENDOSCOPE BRAND	ENDOSCOPE MODEL	COVER SHEATH CODE	
OLYMPUS	ENF-XP	GNF20X30	
	ENF-P4	GNF42X300	
	ENF-GP	GNF42X300	
	ENF-GP2	GNF33X300	
	ENF-V	GNF44X37	
	ENF-V2	GNF42X300	
	ENF-V3	GNF30X30P	
	ENF-V4	GNF30X32O	
	ENF-VH	GNF4230L3	
	ENF-VH2	GNF4230L3	
	ENF-VQ	GNF4230L3	
	STORZ	11101SK-SK2	GNF30X27
11101RP1-2		GNF39X30S	
11101CM		GNF39X30	
11101HD		GNF39X30/39X30V	
11101VP		GNF39X30V	
11102CM		GNF30X30/30X30L	
PENTAX		FNL-7RP3	GNF30X30S
PENTAX	FNL-10RP3	GNF42X300	
	FNL-10RBS	GNF42X300	
	VNL1170K	GNF4430L2	
	VNL1190K	GNF4430L2	
	XION	EF-N SLIM	GNF30X32
		EF-N / EF-NS	GNF39X32
EV-N		GNF44X32	
EV-NC		GNL39X30S	
EV-NE		GNF4430X2	
XN-HD	GNF4232L3		
FENTEX	700710FX	GNF35X30	
MACHIDA	ENT-30PIII	GNF35X30M	
OTOPRONT	VN-S	GNF39X28	
ORL VISION	RS1-RX1	GNFORLV01	
OPTIM	Entity	GNF39X30E	
HENKE SASS WOLF	8300019728	GNF35X30	

C.so Re Umberto, 8 , 10121 Torino (Italy)

Tel **+39 011 0438523**

email:[info@yushuhealthcare.com](mailto:info@yushuhealthcare.com)

[www.yushuhealthcare.com](http://www.yushuhealthcare.com)

