

SCREEN VISION Tissue Transfer & X-Ray System

- Safely image and transport of surgical breast specimens without losing orientation
- Extremely radiolucent and radiopaque grid with pinpointing side indicators
- Precise anatomical assignment of sample margins
- Special stands allowing horizontal and vertical examination

- Once pinned the tumour tissue remains stable and can be easily transported
- On request SCREEN VISION SUPPORT for vertical radiographs




SCREEN VISION

PRODUCT CODE:	SIZE (THICKNESS)
VVX2	20mm
VVX3	30mm
VVX3P	30mm + Pins

(Box of x5 pieces)

Needles with different diameters and lengths may be provided on request.

AREOLAR ESTERNAL DEEP MEDIAL CAUDAL CRANIAL



VVX3P Orientation Pins Version