

RELOCK Preoperative Repositionable Breast Localisation Device

- Preoperative marking of local lesions
- Retractable hook with unique flexibility
- Indicator shows direction of hook
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimetre depth markings assist in needle placement
- Wire fixation clamp to avoid migration
- 9mm loop hook




RELOCK SOFT

PRODUCT CODE:	SIZE
● VVR1905	19G x 5cm
● VVR1907	19G x 7cm
● VVR1910	19G x 10cm
● VVR2007	20G x 7cm
● VVR2010	20G x 10cm
● VVR2015	20G x 15cm

(Box of x10 pieces)



Mammography



Wire Fixation Clamp

Needles with different diameters and lengths may be provided on request.