

INTRO/INTRO PLUS Standard Coaxial Introducer Needle used for Single Access in multiple Biopsies

- Echogenic treatment for ultrasound guided procedures
- Centimetre marking and optional sliding stopper
- Sterilised by ETO, shelf life 5 years
- INTRO (code VVT) with blunt cannula and trocar tip
- INTRO PLUS (code WTC) with sharp cutting cannula and trocar tip



INTRO/INTRO PLUS

PRODUCT CODES:		SIZE
● VVT1105	VVTC 1105	11G x 5cm
● VVT1305	VVTC 1305	13G x 5cm
● VVT1310	VVTC 1310	13G x 10cm
● VVT1315	VVTC 1315	13G x 15cm
● VVT1505	VVTC 1505	15G x 5cm
● VVT1510	VVTC 1510	15G x 10cm
● VVT1515	VVTC 1515	15G x 10cm
● VVT1705	VVTC 1705	17G x 5cm
● VVT1710	VVTC 1710	17G x 10cm
● VVT1715	VVTC 1715	17G x 15cm

(Box of x20 pieces)

Depth marker with code colour



Needles with different diameters and lengths may be provided on request.