


RELOCK DUO Preoperative Repositionable Breast Localisation Device


- Preoperative marking of local lesions
- Retractable hook with unique flexibility
- Echogenic needle tip for ultrasound use
- Hook release at the needle tip for accurate targeting
- Centimetre depth markings assist in needle placement
- Fixation piece for preoperative fixing of the wire
- Twisted wire for enhanced visibility and lower risk of accidental cutting
- Double 5mm loop hook



Half Exposure




Full Exposure




RELOCK DUO

PRODUCT CODE:	SIZE
VVNSP1910	19G x 10cm
VVNSP2005	20G x 5cm
VVNSP2007	20G x 7cm
VVNSP2010	20G x 10cm
VVNSP2012	20G x 12cm
VVNSP2015	20G x 15cm

(Box of x10 pieces)



Mammography



Wire Fixation Clamp

Needles with different diameters and lengths may be provided on request.