

RELOCK PREMIUM DUO Preoperative Repositionable Breast Localisation Device

- Preoperative marking of local lesions
- One-handed design for a convenient manipulation of the ultrasound transducer
- Removable handle for mammography procedure
- Retractable hook with unique flexibility
- Indicator shows direction of hook
- Echogenic needle tip for ultrasound use
- Hooks release at the needle tip for accurate targeting
- Centimetre depth markings assist in needle placement
- Patented stopper for partial hook release prior to confirming position
- Double 5mm loop hook




Half Exposure



Full Exposure



Removable Handle




RELOCK PREMIUM DUO

PRODUCT CODE:	SIZE
VVNSPX1807	18G x 7cm
VVNSPX1810	18G x 10cm
VVNSPX2005	19.5G x 5cm
VVNSPX2009	19.5G x 9cm
VVNSPX2012	19.5G x 12cm

(Box of x10 pieces)



Mammography



Wire Fixation Clamp

Needles with different diameters and lengths may be provided on request.