

GTA (UK) LITEBAG

Drainage bags for collection of fluid during nephrostomy, biliary and abscess drainage.

- **Convenient 600ml size**
 - Not heavy when full
 - Discreet when worn
- **Flock Backing**
 - Eliminates PVC allergy problem
 - High reduction of sweating
 - Optimum patient comfort
- **Pre-attached 50cm (20") silicon tube**
 - Safe and secure connection to drainage catheter
 - Can be cut to desired length
 - Kink resistant
- **Supplied with elastic loop**
 - It can be dressed on the leg, shoulder or waist
- **Anti-reflux valves to avoid retrograde contamination**
- **T-Valve**

GTA (UK) LITEBAG

| PRODUCT CODE: | DESCRIPTION |
|---------------|--------------|
| 500/0601 | T-Valve Cock |

(Box of x25 pieces)



Velcro



T-Valve Cock



The GTA (UK) LiteBag
Is now available on
prescription.